

United States Department of Agriculture
Center for Veterinary Biologics
Testing Protocol

SAM 316

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

Date: November 14, 2014

Number: SAM 316.05

Supersedes: SAM 316.04, June 30, 2011

Standard Requirement: 9 CFR 113.317

Contact: Alethea M. Fry, (515) 337-7200
Peg A. Patterson

Approvals: /s/Geetha B. Srinivas Date: 23Dec14
Geetha B. Srinivas, Section Leader
Virology

/s/Byron E. Rippke Date: 29Dec14
Byron E. Rippke, Director
Policy, Evaluation, and Licensing
Center for Veterinary Biologics

/s/Rebecca L.W. Hyde Date: 29Dec14
Rebecca L.W. Hyde, Section Leader
Quality Management
Center for Veterinary Biologics

United States Department of Agriculture
Animal and Plant Health Inspection Service
P. O. Box 844
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

Table of Contents

- 1. Introduction**
- 2. Materials**
 - 2.1 Equipment/instrumentation**
 - 2.2 Reagents/supplies**
- 3. Preparation for the Test**
 - 3.1 Personnel qualifications/training**
 - 3.2 Preparation of equipment/instrumentation**
 - 3.3 Preparation of reagents/control procedures**
 - 3.4 Preparation of the sample**
- 4. Performance of the Test**
- 5. Interpretation of the Test Results**
- 6. Report of Test Results**
- 7. References**
- 8. Summary of Revisions**

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

1. Introduction

This Supplemental Assay Method (SAM) describes an *in vitro* test method for assaying modified-live canine parvovirus (CPV) vaccines for viral content. The method uses the Crandell feline kidney (CRFK) cell line as the test system. CPV endpoint is determined by staining inoculated cell cultures by an indirect fluorescent antibody (IFA) staining method.

2. Materials

2.1 Equipment/instrumentation

Equivalent equipment or instrumentation may be substituted for any brand name listed below.

2.1.1 Incubator, $36^{\circ} \pm 2^{\circ}\text{C}$, high humidity, $5\% \pm 1\%$ CO_2 (Model 3336, Forma Scientific Inc.)

2.1.2 Incubator aerobic (Model 2, Precision Scientific)

2.1.3 Water bath

2.1.4 Microscope, ultraviolet (UV) light (Model BH2, Olympus America)

2.1.5 Vortex mixer (Vortex-2 Genie, Model G-560, Scientific Industries)

2.1.6 Micropipettor, tips, and/or motorized microliter pipette or equivalent and tips

2.1.7 Microscope slide glass staining dish with rack (glass staining dish)

2.2 Reagents/supplies

Equivalent reagents or supplies may be substituted for any brand name listed below. All reagents and supplies must be sterile.

2.2.1 CPV Positive Control, KB5 strain

2.2.2 CRFK cell culture, free of extraneous agents as tested by title 9, *Code of Federal Regulations* (9 CFR)

2.2.3 Feline panleukopenia virus antiserum (FPV Antiserum) or Canine Parvovirus antiserum (CPV Antiserum)

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

2.2.4 Minimum essential medium (MEM)

1. 9.61 g MEM with Earles salts without bicarbonate
2. 1.1 g sodium bicarbonate (NaHCO_3)
3. Dissolve with 900 mL deionized water (DI).
4. Add 5.0 g lactalbumin hydrolysate or edamine to 10 mL DI. Heat to $60^\circ \pm 2^\circ\text{C}$ until dissolved. Add to **Step 3** with constant mixing.
5. Q.S. to 1000 mL with DI; adjust pH to 6.8-6.9 with 2N hydrochloric acid (HCl).
6. Sterilize through a 0.22- μm filter.
7. Aseptically add 50 $\mu\text{g}/\text{mL}$ gentamicin sulfate
8. Store at $2^\circ - 7^\circ\text{C}$.

2.2.5 Growth Medium

1. 920 mL of MEM
2. Aseptically add:
 - a. 70 mL gamma-irradiated fetal bovine serum (FBS)
 - b. 10 mL L-glutamine (200 mM)

2.2.6 Dulbecco's phosphate buffered saline (DPBS)

1. 8.0 g sodium chloride (NaCl)
2. 0.2 g potassium chloride (KCl)
3. 0.2 g potassium phosphate, monobasic, anhydrous (KH_2PO_4)
4. 0.1 g magnesium chloride, hexahydrate ($\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$)
5. Dissolve reagents with 900 mL DI.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

6. Dissolve 1.03 g sodium phosphate, dibasic, anhydrous (Na_2HPO_4) with 10 mL DI, heat to $60^\circ \pm 2^\circ\text{C}$ until dissolved; add to **Step 5** with constant mixing.

7. Dissolve 0.1 g calcium chloride, anhydrous (CaCl_2) with 10 mL DI; add slowly to **Step 6** to avoid precipitation.

8. Q.S. to 1000 mL with DI, adjust pH to 7.0-7.3 with 2N HCl.

9. Sterilize through a 0.22- μm filter.

2.2.7 Cell culture glass slides, 8-well (Lab-Tek[®] Slides)

2.2.8 Polystyrene tubes, 12 x 75-mm

2.2.9 Appropriate anti-species IgG (H&L) fluorescein isothiocyanate labeled conjugate (Anti-species Conjugate)

2.2.10 100% acetone

2.2.11 Syringe, 1-mL and needle, 20-gauge x 1 1/2-inch

2.2.12 Pipette-Aid

2.2.13 Disposal transfer pipette, 4.6-mL

3. Preparation for the Test

3.1 Personnel qualifications/training

Personnel shall have experience in the preparation and maintenance of cell culture as well as in the propagation of animal viruses and the quantitation of virus infectivity by the IFA staining method.

3.2 Preparation of equipment/instrumentation

3.2.1 On the day of inoculation, set a water bath at $36^\circ \pm 2^\circ\text{C}$.

3.2.2 On the day of the IFA test, prepare a humidity chamber in the aerobic incubator by filling a pan in the bottom with DI.

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

3.3 Preparation of reagents/control procedures

3.3.1 Preparation of CRFK Slides

1. Cells are prepared from healthy, confluent CRFK cells, that are maintained by passing every 3 to 4 days. On the day of test initiation, add 0.4 mL/well of $10^{5.2}$ to $10^{5.4}$ cells/mL diluted in Growth Media into all wells of a Lab-Tek® slide. Prepare sufficient Lab-Tek® Slides to allow 25 wells for controls and 20 wells for each Test Vaccine. Incubate at $36^{\circ}\pm 2^{\circ}\text{C}$ in a CO_2 incubator. These become the CRFK Slides.

2. Use seeded CRFK Slides within 4 hours.

3.3.2 Preparation of the CPV Working Positive Control

1. On the day of inoculation, rapidly thaw a vial of CPV Positive Control in a $36^{\circ}\pm 2^{\circ}\text{C}$ water bath.

2. Dispense 1.8 mL of MEM into each of 7, 12 x 75-mm polystyrene tubes labeled 10^{-1} through 10^{-7} .

3. Transfer 200 μL of the CPV Positive Control to the tube labeled 10^{-1} ; discard pipette tip. Mix by vortexing.

4. Using a new pipette tip, transfer 200 μL from the 10^{-1} labeled tube (**Step 3**) to the 10^{-2} tube; discard pipette tip. Mix by vortexing.

5. Repeat **Step 4** for each subsequent dilution, transferring 200 μL from the previous dilution to the next dilution tube until the tenfold dilution series is completed.

3.3.3 Preparation of Working FPV or CPV Antiserum

On the day of the IFA test, dilute FPV or CPV Antiserum in DPBS to the IFA working dilution on the CVB Reagent Data Sheet or as determined for that specific antiserum.

3.3.4 Preparation of Working Anti-species Conjugate

On the day of the IFA test, dilute Anti-species Conjugate in DPBS to the working dilution according to the manufacturer's recommendations.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

3.4 Preparation of the sample

3.4.1 The initial test of a Test Vaccine will be with a single vial (a single sample from 1 vial). On the day of inoculation, rehydrate a vial of the Test Vaccine by transferring 1.0 mL for a 1-mL-dose vaccine, 0.5 mL for 1/2-mL-dose vaccines, etc., of the provided diluent into the vial containing the lyophilized Test Vaccine. Use a sterile 1.0-mL syringe and an 18-gauge x 1 1/2-inch needle; mix by vortexing. Incubate for 15 ± 5 minutes at room temperature.

3.4.2 Dispense 1.8 ml MEM into each of 6, 12 x 75-mm polystyrene tubes labeled 10^{-1} through 10^{-6} .

3.4.3 With a micropipettor, transfer 200 μ L from the rehydrated Test Vaccine vial to the tube labeled 10^{-1} ; discard the pipette tip. Mix by vortexing.

3.4.4 Transfer 200 μ L from the 10^{-1} labeled tube (**Section 3.4.3**) to the 10^{-2} tube; discard the pipette tip. Mix by vortexing.

3.4.5 Repeat **Section 3.4.4** for each subsequent dilution, transferring 200 μ L from the previous dilution to the next dilution tube until the tenfold dilution series is completed.

4. Performance of the Test

4.1 Inoculate 5 wells/dilution of a CRFK Slide with 100 μ L/well from dilutions 10^{-6} through 10^{-3} of the Test Vaccine. Tip changes are not necessary between each dilution in a series if pipetting from the most dilute to the most concentrated within that series (e.g., 10^{-6} through 10^{-3}).

4.2 Inoculate 5 wells/dilution of a CRFK Slide with 100 μ L/well, from dilutions 10^{-6} through 10^{-3} of the CPV Working Positive Control. Tip changes are not necessary between each dilution in a series if pipetting from the most dilute to the most concentrated within that series (e.g., 10^{-6} through 10^{-3}).

4.3 Five uninoculated wells serve as a Negative Cell Control.

4.4 Incubate CRFK Slides in a $36^{\circ} \pm 2^{\circ}\text{C}$ CO_2 incubator for 5 days \pm 1 day.

4.5 Following incubation, decant the media from the CRFK Slide and remove the plastic walls by twisting them away from the CRFK Slide, leaving the gasket attached to the slide.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

- 4.6** Place the CRFK Slides in a slide rack; place the rack in a glass staining dish filled with DPBS. Let stand for 15 ± 5 minutes at room temperature.
- 4.7** Discard the DPBS; fix the CRFK Slides in 100% acetone for 15 ± 5 minutes at room temperature. Remove and allow to air dry.
- 4.8** Pipette 75 ± 25 μ L of the Working FPV Antiserum into each well of the CRFK Slides. Incubate for 30 ± 5 minutes in the aerobic incubator at $36^{\circ} \pm 2^{\circ}$ C.
- 4.9** Wash per **Section 4.6**.
- 4.10** Pipette 75 ± 25 μ L of the Working Anti-species Conjugate into each well of the CRFK Slides. Incubate for 30 ± 5 minutes in the aerobic incubator at $36^{\circ} \pm 2^{\circ}$ C.
- 4.11** Wash per **Section 4.6**. Discard the DPBS.
- 4.12** Rinse the CRFK Slides with DI; allow to air dry.
- 4.13** Read at 100-200X magnification with a UV-light microscope; examine the cell monolayer for typical CPV apple-green nuclear fluorescence.
- 4.13.1** Wells containing 1 or more cells with specific CPV fluorescence are considered to be positive.
- 4.13.2** Results are recorded as number of IFA positive wells versus total number of wells examined for each dilution of a Test Vaccine and the CPV Working Positive Control.
- 4.14** Calculate the 50% fluorescent antibody infective dose (FAID₅₀) of the Test Vaccine and the CPV Working Positive Control using the method of Spearman-Kärber as modified by Finney.

Example:

10^{-3} dilution of Test Vaccine = 5/5 wells IFA positive
 10^{-4} dilution of Test Vaccine = 4/5 wells IFA positive
 10^{-5} dilution of Test Vaccine = 2/5 wells IFA positive
 10^{-6} dilution of Test Vaccine = 0/5 wells IFA positive

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

Spearman-Kärber formula:

Test Vaccine Titer = $(X - d/2 + [d \cdot S])$ where:

X = \log_{10} of highest dilution with all wells IFA positive (3)

d = \log_{10} of tenfold dilution factor (1)

S = sum of proportions of wells IFA positive for all dilutions tested:

$$\frac{5}{5} + \frac{4}{5} + \frac{2}{5} + \frac{0}{5} + \frac{11}{5} = 2.2$$

Test Vaccine titer = $(3 - 1/2) + (1 \cdot 2.2) = 4.7$

Adjust the titer to the Test Vaccine dose as follows:

A. divide the **Test Vaccine Dose** by the **Inoculation Dose**

Test Vaccine Dose = manufacturer's recommended vaccination dose (for this CPV test vaccine, the recommended dose is 1 mL)

Inoculation Dose = amount of diluted Test Vaccine added to each well of the Test Slide (for this CPV test vaccine, the inoculation dose is 0.1 mL)

$$\frac{1 \text{ mL dose}}{0.1 \text{ mL inoculum}} = 10$$

B. calculate \log_{10} of value in **A** and add it to the **Test Vaccine titer** as illustrated below:

$$\text{Log of } 10 = 1.0$$

Test Vaccine titer = $4.7 + 1.0 = 5.7$

Therefore the titer of the **CPV Test Vaccine** is $10^{5.7}$ FAID₅₀/mL.

5. Interpretation of the Test Results

5.1 Validity requirements:

5.1.1 The calculated FAID₅₀ titer of the CPV Working Positive Control must fall within plus or minus 2 standard deviations (± 2 SD) of its mean titer, as established from a minimum of 10 previously determined titers.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

5.1.2 The lowest inoculated dilution of the CPV Working Positive Control must exhibit a 100% positive IFA reaction (5/5). If an endpoint is not reached (1 or more wells are IFA positive at the highest dilution), the titer is expressed as “greater than or equal to” the calculated titer. If an endpoint is critical to testing, the highest (most dilute) must exhibit no positive IFA reaction (0/5).

5.1.3 The Uninoculated Cell Control must not exhibit any cytopathic effect, specific CPV fluorescence, or cloudy media that would indicate contamination.

5.2 If the validity requirements are not met, then the assay is considered a **NO TEST** and may be retested without prejudice.

5.3 If the validity requirements are met and the titer of the Test Vaccine is greater than or equal to the titer contained in the Animal and Plant Health Inspection Service (APHIS) filed Outline of Production for the product under test, the Test Vaccine is considered **SATISFACTORY**.

5.4 If the validity requirements are met, but the titer of the Test Vaccine is less than the required minimum titer contained in the APHIS filed Outline of Production for the product under test, the Test Vaccine may be retested in accordance with the 9 CFR 113.8.

6. Report of Test Results

Report results as the FAID₅₀ per dose of the Test Vaccine.

7. References

7.1 Title 9, *Code of Federal Regulations*, part 113.317, U.S. Government Printing Office, Washington, DC.

7.2 Cottral, GE, (Ed.), 1978. *Manual of standardized methods for veterinary microbiology*. Comstock Publishing Associates, Ithaca, New York. pg. 731.

7.3 Finney, DJ, 1978. *Statistical methods in biological assay*. Griffin, London. 3rd edition, pp. 394-401.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

8. Summary of Revisions

Version .05

- The Contact information has been updated; however, the Virology Section has elected to keep the same next review date for the document.

Version .04

- The Contact information has been updated.

Version .03

- The term “Reference” has been changed to “Positive Control” throughout the document.

Version .02

This document was revised to clarify the practices currently in use in the Center for Veterinary Biologics and to provide additional detail. While no significant changes impacting the test were made from the previous protocol, the following changes have been made:

- **2.1.2, 2.1.3** The temperature setting has been removed.
- **2.1.6, 2.2.12, 2.2.13** Equipment has been added.
- **2.2.3, 3.3.3** The reagent name has been corrected to “or CPV antiserum”.
- **2.2.4.2** The amount of sodium bicarbonate (NaHCO₃) has been changed from 2.2 g to 1.1 g.
- **2.2.4.7** Penicillin and streptomycin have been deleted.
- **2.2.9, 3.3.4, 4.10** The conjugate name has been corrected to Anti-Species IgG.
- **3.1** Personnel experience has been clarified.
- **3.3.1** The cell count has been clarified.
- **4.2** Inoculation of the Reference has been clarified.
- **4.4** The incubation has been changed from hours to days.
- **4.9, 4.11** The wash step has been corrected.

UNCONTROLLED COPY

Supplemental Assay Method for the Titration of Canine Parvovirus in Cell Culture

- **4.13.2** The recording step has been corrected.
- **4.14** Additional steps have been added to clarify the titer calculations by the Spearman-Kärber formula.
- **5.1.2** Recording the rate of positive reaction for validity requirements.
- “Reference and Reagent Sheet” has been changed to “Reagent Data Sheet” throughout the document.
- The refrigeration temperatures have been changed from $4^{\circ}\pm 2^{\circ}\text{C}$ to $2^{\circ}\text{-}7^{\circ}\text{C}$. This reflects the parameters established and monitored by the Rees system.
- The word “chamber” has been changed to “well” throughout the document.
- The footnotes have been deleted with any pertinent references now noted next to the individual items.
- “Test Serial” has been changed to “Test Vaccine” throughout the document.

UNCONTROLLED COPY