

Diagnosics of *Dickey dianthicola*

These protocols are for detection and identification of *Dickeya* sp., particularly *D. dianthicola*, in plant samples. Please note that these methods have not been validated by PPQ. They are the results of a literature survey and experience.

Isolation from tissue

Sampling from potato tuber

1. Gently wash tubers with 10% bleach to remove excess soil. Rinse with water and damp dry.
2. Using a disposable razor blade, cut into the stolon tip of the tuber and look for infected tissue. Take sample tissue from the margins of the infected area (see image below). Including both necrotic and healthy tissue in the sample, cut enough material for approximately 1.0 g of tissue. Chop into pieces on sterile weigh boat.
3. Place ~100 mg of the sample into a Q-Biogene Lysing matrix tube to be used for nucleic acid extraction using the Qiagen DNeasy Plant Mini Kit.
4. Place ~500 mg into a mesh extraction bag (e.g., Agdia or Bioreba) with 3 mL of 50 mM phosphate buffer (PB, pH 7.2, containing 4.26 g Na₂HPO₄ and 2.72 g KH₂PO₄/liter) and thoroughly crush sample.
 - a. Pipette out 500 µL of liquid and place into 1.5 mL microcentrifuge tube to be used for bacterial isolation.

Sampling from potato stem

1. Wipe external part of stem with 10% bleach wipe to remove excess soil. Rinse with water and damp dry.
2. Using a disposable razor blade, remove a ~5 cm section of stem containing both diseased and healthy tissue, near the crown region if possible.
3. Chop the tissue using a new disposable razor blade on a sterile weight boat.
4. Place ~100 mg of tissue into a Q-Biogene Lysing matrix to be used for nucleic acid extraction using the Qiagen DNeasy Plant Mini Kit.
5. Place ~500 mg of the tissue into a mesh extraction bag with 3 ml 50 mM PO₄ buffer and crush sample.
 - a. Pipette out 500 µL of liquid and place into 1.5 mL microcentrifuge tube to be used for bacterial isolation.

DNA Extraction

Thoroughly disrupt tissue in a lysing matrix tube using a FastPrep machine, bead-beater, or similar. Extract DNA according to Qiagen instructions.




Figure 1: Locations circled for infected tuber and stem sampling.



Bacterial Isolation

1. Allow crushed material to soak in buffer for 5-10 minutes at room temperature.
2. Using a sterile inoculating loop, take one loop of suspension and streak onto three petri plates:
 - a. Tryptic Soy Agar (TSA)
 - b. Nutrient Agar (NA)
 - c. Crystal violet pectate medium (CVP) [Sonia N. Humphris et al., *Diagnosis of *Pecobacterium* and *Dickeya* spp.* 2015]
3. Optional: Pipette 20-40 μL of the sample suspension onto three fresh petri plates (TSA, NA, and CVP) and streak the liquid fully across the media surface using a sterile plastic spreader. Store sample suspension at 4°C.
4. Seal each plate with parafilm. Incubate at 32°-37°C for 24 hours. Store sample suspension at 4°C.
5. Examine plates for single colonies, ideally about 2 mm or larger in diameter.
6. Select colony using a 1 μL sterile plastic loop and place into 1.5 mL microcentrifuge tube with 200 μL 10% DMSO 100 mM Tris.
7. Incubate samples at 95°C for 2 min. Immediately place on ice.
8. Store at 4°C for short term, and -20°C for long term storage.

Figure 2: *Dickeya* colonies on select media for bacterial isolation. **A)** CVP **B)** TSA



Conventional PCR Assays

Master Mixes for PCR

1. Keep all reagents and samples on ice when working.
2. Prepare Master Mix following Table 1 and aliquot 23 μL into each 0.2 mL PCR tube.
3. Add 2.0 μL of undiluted sample DNA for a total reaction volume of 25 μL . Each PCR run should include two controls: *Dickeya sp.* control and non-template control (NTC-MG water)

Table 1: For all PCR reactions use the following master mix:

Reagent	Amount for 25 μL reaction volume
Nuclease-Free H ₂ O	9.5 μL
GoTaq® Green Master Mix (Promega M712)	12.5 μL
6 μM Primer mix (Forward+Reverse)	1.0 μL

Diagnostics

NOTE: although contemporary publications describing bacterial characterization using MLSA (multi-locus sequence analysis) describe using DNA from pure bacterial colonies, we have had good success targeting the genetic locations in the following table (*fliC*, *dnaX*, and *gyrB*) using total DNA extracted using a Qiagen kit which has tested positive for presence of *Dickeya* using one or both of the PCR screening tests (DF/DR, *pelADE1/pelADE2*). When successful, this method vastly reduces the time and labor required to determine a sample positive for *Dickeya*. Occasionally, more than one bacterial species will be amplified, which results in an unreadable sequence when direct sequencing is used. For this reason, we recommend testing with all three MLSA primer sets in order to obtain at least 2 sequences with the required homology to *D. dianthicola*.

1. To screen suspect samples for *Dickeya sp.*, run primer sets *pelADE1/2* and/or DF/DR conventional PCR on all samples.
2. Using samples that gave positive results for either or both *pelADE1/2* / DF/DR tests, continue testing using conventional PCR primer sets *flic1/2*, *dnaXF/R*, and *gyrBF/R*.
3. For each run, use 10 μL of PCR product on 1.5% agarose gel using 1X TAE buffer. Run at 100V for 60 min. and stain with ethidium bromide for visualization.

Table 2: Conventional PCR primers

Primer Name	Sequence 5'-3'	Amplicon Size (bp)
<i>pelADE1</i>	GAT CAG AAA GCC CGC AGC CAG AT	420
<i>pelADE2</i>	CTG TGG CCG ATC AGG ATG GTT TTG TCG TGC	
DF	AGA GTC AAA AGC GTC TTG	133
DR	TTT CAC CCA CCG TCA GTC	
<i>flic1</i>	TAT CAA CAG CGC CAA AGA CAA CGC	650
<i>flic2</i>	ACG GCT CAT GTT GGA TAC TTC GTT	
<i>dnaXF</i>	TAT CAG GTY CTT GCC CGT AAG TGG	535
<i>dnaXR</i>	TCG ACA TCC ARC GCY TTG AGA TG	
<i>gyrBF</i>	TAA GTT YGA CGA YAA CTC STA YAA RGT	974
<i>gyrBR</i>	CCC CTT CCA CCA GGT ASA GTT C	

Table 3: Conventional PCR cycling conditions

Primers	Step 1	Step 2	Step 3
pelADE1/2	94°C for 2 min	34 cycles: 94°C for 30s, 72 °C for 90s	72 °C for 5 min
DF/DR	94°C for 2 min	30 cycles: 94°C for 20s, 62°C for 90s, 72 °C for 90s	72 °C for 5 min
fliC1/2	95°C for 1 min	35 cycles: 95°C for 20s, 55°C for 30s, 72 °C for 45s	72 °C for 5 min
dnaXF/R	94°C for 3 min	35 cycles: 94°C for 60s, 59°C for 60s, 72 °C for 120s	72 °C for 5 min
gyrBF/R	94°C for 4 min	35 cycles: 94°C for 60s, 56°C for 60s, 72 °C for 120s	72 °C for 10 min

Evaluation of Samples




Figure 3: pelADE1/pelADE2: Conventional PCR assay of *Dickeya* suspect samples. Lanes 2, 10 and 11 are positive.




Figure 4: DF/DR: Conventional PCR assay of suspect *Dickeya* samples. Lanes 1, 2, 3, 5, and 6 are positive.




Figure 5: dnaXF/dnaXR: Conventional PCR assay of suspect *Dickeya* samples. Lanes 2, 3, and 4 are positive.




Figure 6: gyrBF/gyrBR: Conventional PCR assay of suspect *Dickeya* samples. Lanes 1, 2, and 3 are positive.




Figure 7: *fliC1/fliC2*: Conventional PCR assay of suspect *Dickeya* samples. Lanes 1, 2, and 3 are positive.

Sequence analysis of directly sequenced PCR products

1. PCR products from *fliC*, *dnaX*, and *gyrB* that show a positive result matching the positive control with no contamination in the NTC can be directly sent for sequencing. Please note that presence of a PCR product of the apparent correct size is not proof that *Dickeya* DNA has been amplified.
2. In order to determine positive for *Dickeya*, BLAST aligned consensus sequences against the following accession # and % match.
 - a. *fliC*: $\geq 98\%$ homology with accession KF442436.1: positive for *D. dianthicola*
 - b. *dnaX*: $\geq 98\%$ homology with accession LC043164.1: positive for *D. dianthicola*
 - c. *gyrB*: $\geq 98\%$ homology with accession AOOK01000040.1: positive for *D. dianthicola*

References:

Characterization of *Erwinia chrysanthemi* by Pectinolytic Isozyme Polymorphism and Restriction Fragment Length Polymorphism Analysis of PCR-Amplified Fragments of *pel* Genes. 1996. Nassar, A., Darrasse, A., Lemattre, M., Kotoujansky, A., Dervin, C., Vedel, R., and Bertheau, Y. 1996. *Applied and Environmental Microbiology* 62: 2228–2235.

Biochemical and genetical analysis reveal a new clade of biovar 3 *Dickeya* spp. strains isolated from potato in Europe. 2009. Slawiak, M., van Beckhoven, J., Speksnijder, A., Czajkowski, R., Grabe, G., and van der Wolf, J. *Eur J Plant Pathol* 125:245–261.

Sequence Diversity in the *Dickeya fliC* Gene: Phylogeny of the *Dickeya* Genus and TaqMan PCR for *D. solani*, New Biovar 3 Variant on Potato in Europe. 2012. Van Vaerenberg, J., Baeyen, S., De Vos, P., and Maes, M. *PLoS ONE* 7, e35738.

Symptoms and yield reduction caused by *Dickeya* spp. strains isolated from potato and river water in Finland. 2009. Laruila, J., Hannukkala, A., Nykyri, J., Pasanen, M., Helias, V., Garland L., Pirhonen, M. Eur J Plant Pathol 126:249-262.

Phylogeny and classification of *Dickeya* based on multilocus sequence analysis. 2013. Marrero, G., Schneider, KL., Jenkins, DM., Alvarez, AM. Int J Syst Evol Microbiol (pt 9):3524-39.

Detection of the Bacterial Potato Pathogens *Pectobacterium* and *Dickeya* spp. Using Conventional and Real-Time PCR. 2015. Humphris, S., Cahill, G., Elphinstone, J., Kelly, R., Parkinson, N., Pritchard, L., Toh, I., Sandler, G. Plant Pathology Techniques and Protocols, Second edition . Vol 1302.