

Animal and Plant Health Inspection Service

Veterinary Services

NAHMS Goat 2019 Parasite Pre-deworming (Kit A) Collection Record

National Animal Health Monitoring System

2150 Centre Ave, Bldg B Fort Collins, CO 80526

Form Approved OMB Number 0579-0354 Approval expires:04/30/2022

Kit contents:

25 small Whirl-Pak® bags, lubricant, 2 ice packs, 1 liner bag, 1 medium insulated cooler, and paperwork that includes submission form, labels, and 1 UPS airbill addressed to LSU in Baton Rouge, LA. You will need to provide you own gloves. Use clean gloves for each goat.

Collection Instructions

- 1. Collect samples Sunday-Wednesday. Sample goats that have not been dewormed in the previous 60 days. We recommend deworming animals at the time of collection.
- 2. The number of samples collected is based on the number of resident goats on the operation. Goats sampled should represent the goats and kids on the operation in terms of age, sex, breed, and use. We recommend including goats that the owner believes are likely to have worms. Use the following chart for determining sample numbers:

Number of Goats on an Operation	Sample Number
1 to 19 goats	Sample all goats
20 to 49 goats	Sample 20 goats
50 or more goats	Sample 25 goats

- 3. From each goat, collect 5-6 fecal pellets. Collect samples from the rectum when possible. Rectal retrieval might not be possible on some goats (e.g. preweaned kids), and fresh off the ground samples are acceptable. On each label, write the goat's name or ID and attach the label to the sample bag.
- 4. Samples must be fresh (not petrified). Do not exclude diarrhea samples.
- 5. If the sample cannot be associated with a specific goat, write NO INDIV GOAT in column for name on ID, but complete as much of the other information as possible. However, only identified goats with a pre-deworming sample submission will be tested post-deworming.
- 6. If deworming at time of collection, please include the used dewormer tube, label, or insert in the sample shipping box that is sent to the lab.
- 7. Cool samples down as soon as possible (in a refrigerator or cooler).

According to the Paperwork Reduction Act of 1995, an agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a valid OMB control number. The valid OMB control number for this information collection is 0579-0354. The time required to complete this information collection is estimated to average 2.0 hours per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information.

NAHMS-456 JUN 2019

RECTAL RETRIEVAL

To avoid contamination from common organisms on the ground, rectal retrieval is best. Rectal retrieval might not be possible on some goats (e.g. preweaned kids), and fresh off the ground samples are acceptable.



- 1. Apply lubricating jelly to the glove before entering the rectum.
 - Lightly stroking the rectum might encourage defecation.



2. Retrieve a minimum of 5-6 pellets per animal.

- 3. Place pellets in a Whirl-Pak® bag, removing any excess air. On each label, write the goat's name or ID and attach it on to the bag.
- 4. Continue collecting samples from other goats using a clean glove for each.

Collection Record Form Instructions

- 1. **Using a ballpoint pen,** record samples on the appropriate lines and complete all information requested.
- 2. Send the original white and yellow collection records to the lab and leave the pink copy with the producer. Place the producer's pink collection record copy in the post-deworming kit so that the producer can reference this copy and match IDs when they collect the post-deworming samples.

Shipping Instructions

- 1. **Ship on Monday-Wednesday.** Keep samples cool and ship within 24 hours of collection. Wednesday collections must be shipped the same day. Do not collect or ship samples Thursday through Saturday.
- 2. Place all the samples in the liner bag and tie shut. Place an ice pack on the top and bottom of the samples. Add filler to box if necessary. Close the insulated cooler box and place the white and yellow collection records on top of the cooler box lid. The pink copy stays with the producer.
- 3. Secure the box and ship to LSU, in Baton Rouge, LA, within 24 hours. A shipping airbill is provided in the kit. **Ship only Monday-Wednesday.**

NOTE: Remove or black out all extraneous labels on outside of box.

NAHMS ID #:	á	aid

Primary collector name/phone: aname/aph Number of goats on operation: _____anum

Total	sampling	time:		ahr
404		400	400	

Parasite Pre-deworming Kit A Collection Record

a101	a102	a103	a104	a105	a106	a107	a108	a109a-b	a110a-b	a111
Sample #	Goat name or ID	Age (months or years)	Goat Type 1=pregnant doe 2=nursing doe 3=preweaned kid 4=weaned kid 5=open doe 6=buck 7=wether	If goat type= 1 or 2 provide date kidded or expected to kid	Breed [See codes below]	FAMACHA score	# of times dewormed in last 12 months [Not including today]	Dewormer used at last deworming prior to this study [Enter codes from reference card]	Dewormer used at time of this fecal collection [Enter codes from reference card]	Body condition score 1=thin 2=normal 3=fat
1		mo OR yr								
2		mo OR yr								
3		mo OR yr								
4		mo OR yr								
5		mo OR yr								
6		mo OR yr								
7		mo OR yr								
8		mo OR yr								
9		mo OR yr								
10		mo OR yr								

Breed Codes:

1=Alpine 4=Cashmere 2=Angora 5=Fainting goats 3=Boer 6=Kiko 9=Nubian

7=LaMancha 8=Nigerian dwarf

10=Oberhasli 11=Pygmy

12=Pygora

13=Saanen 14=Sable

15=Savannah

16=Spanish 17=Toggenburg

18=Crossbred (specify_____)

19=Other (specify __

Sample #	Goat name or ID	Age (months or years)	Goat Type 1=pregnant doe 2=nursing doe 3=preweaned kid 4=weaned kid 5=open doe 6=buck 7=wether	If goat type= 1 or 2 provide date kidded or expected to kid	Breed [See codes below]	FAMACHA score	# of times dewormed in last 12 months [Not including today]	Dewormer used at last deworming prior to this study [Enter codes from reference card]	Dewormer used at time of this fecal collection [Enter codes from reference card]	Body condition score 1=thin 2=normal 3=fat
11		mo OR yr								
12		mo OR yr								
13		mo OR yr								
14		mo OR yr								
15		mo OR yr								
16		mo OR yr								
17		mo OR yr								
18		mo OR yr								
19		mo OR yr								
20		mo OR yr								

Breed Codes:

1=Alpine4=Cashmere7=LaMancha10=Oberhasli13=Saanen16=Spanish2=Angora5=Fainting goats8=Nigerian dwarf11=Pygmy14=Sable17=Toggenburg

3=Boer 6=Kiko 9=Nubian 12=Pygora 15=Savannah 18=Crossbred (specify______

19=Other (specify _____)

Sample #	Goat name or ID	Age (months or years)	Goat Type 1=pregnant doe 2=nursing doe 3=preweaned kid 4=weaned kid 5=open doe 6=buck 7=wether	If goat type= 1 or 2 provide date kidded or expected to kid	Breed [See codes below]	FAMACHA score	# of times dewormed in last 12 months [Not including today]	Dewormer used at last deworming prior to this study [Enter codes from reference card]	Dewormer used at time of this fecal collection [Enter codes from reference card]	Body condition score 1=thin 2=normal 3=fat
			Continue collec	tion if there are s	0 or more doe	s and kids on	the operation.			
21		mo OR yr								
22		mo OR vr								
23		mo OR yr								
24		mo OR yr								
25		mo OR yr								

Send the original white and yellow collection records to the lab.

Leave the pink copy with the producer.

Breed Codes:

1=Alpine4=Cashmere7=LaMancha10=Oberhasli13=Saanen16=Spanish2=Angora5=Fainting goats8=Nigerian dwarf11=Pygmy14=Sable17=Toggenburg

3=Boer 6=Kiko 9=Nubian 12=Pygora 15=Savannah 18=Crossbred (specify______

19=Other (specify _____)