United States Department of Agriculture Animal and Plant Health Inspection Service Center for Veterinary Biologics P. O. Box 844 Ames, IA 50010

1. Reagent Name: Negative Guinea Pig Serum

2. Strain or Source: Not Applicable

3. Lot Number: IRP 550 (07)

4. Fill Date: April 19, 2007

5. Expiration Date: December 31, 2022

Precautions: There are no known hazards associated with the use of this reagent.

- **6. Intended Use:** IRP 550 serves as the negative control serum when quantitating the tetanus antitoxin titer of guinea pig serum by ELISA.
- 7. Instructions for Use: When testing the potency of guinea pig tetanus antitoxin by ELISA, as described in **Supplemental Assay Method (SAM) 217**, dilute the serum 1:4,000 in PBS/milk/Tween diluent (0.015 M phosphate buffered saline, pH 7.2, 1% nonfat dry milk, 0.05% Tween 20). The test dilution is prepared by adding 0.1 mL of well mixed IRP 550 to 9.9 mL of PBS/milk/Tween diluent. The serum is further diluted to 1:4,000 by adding 0.1 mL of the 1:100 dilution to 3.9 mL of PBS/milk/Tween diluent.

8. Test of Reagent:

Determination of Absorbance Values - Absorbance values of less than 0.05 were recorded when IRP 550 was tested by ELISA as described in **SAM 217**.

Sterility test - The serum was tested for sterility and found to be free of viable bacteria and fungi according to procedures outlined in title 9, *Code of Federal Regulations* (9 CFR), section 113.25.

- **9. Container Size, Type, Weight, or Volume:** Two-mL glass vials containing 1.2 mL of serum.
- **10. Storage Conditions:** Store at -20° C or lower. Once the serum has been thawed store at $4^{\circ} \pm 2^{\circ}$ C.
- **11. CVB Technical Contact:** Bacteriology Section, Center for Veterinary Biologics, (515) 337-6140 or FAX (515) 337-7673.

BBDAT0113.04 Page 1 of 2

- 12. Origin and Passage History: Not applicable.
- **13. Method of Preparation:** Serum was collected from 500-800 gram nonvaccinated female and male guinea pigs (purchased from Small Stock and Charles River Laboratories). The serum was passed through a sterile Millipore filtration unit containing a 0.22-μm membrane, dispensed into sterile vials, and stored at -20°C.

14. Other: None

Reagent orders and feedback should be sent *including phone number* to the following email address: CVB@aphis.usda.gov

Reagent orders forms (APHIS 2018) are available from: https://www.aphis.usda.gov/library/forms/pdf/APHIS_2018.pdf

REVISED: 19Jan18 jmw

BBDAT0113.04 Page 2 of 2