

## Spray-Dried Porcine Plasma to Brazil

This protocol is intended to cover spray-dried porcine plasma intended for use in animal feed exported to Brazil.

### **General Instructions for Exporters**

1. Use the fillable letterhead certificate: “INTERNATIONAL VETERINARY CERTIFICATE FOR THE EXPORT OF SPRAY DRIED PORCINE BLOOD PRODUCTS FOR ANIMAL FEEDING” available on the main Brazil IREG page.
2. Have the manufacturer prepare a notarized affidavit with the required certification statements noted under the notarized affidavit line in the “Declarations” section of the letterhead certificate. If assistance is needed in preparing a proper notarized affidavit, please contact the pertinent Veterinary Services ([VS Service Center](#)) that will be endorsing the export certificate.
3. For Statement #3, select the most appropriate option(s); if more than one treatment applies, select all that apply.
4. For Statement #4, select the appropriate option. Only one option may be selected. Note that for the third option, lot-specific results demonstrating the absence of ruminant protein by PCR must be submitted to VS with the certificate.
5. Complete the “PORT” and “DESCRIPTION OF CONSIGNMENT” section of the letterhead certificate by following the [general instructions](#) for the VS Form 16-4. If assistance is needed, please contact the pertinent VS Service Center that will be endorsing the export certificates. In addition to the basic information required, Brazil also requires the following:
  - a. Product: type of product, product name, species of origin, type and number of packaging, net weight
  - b. Identification: lot number(s), date or range of dates of production, container and seal number (for bulk shipments)
6. Submit the certificate and affidavit, including negative ruminant PCR results if required, to the pertinent VS Service Center
7. The VS Service Center will print the certificate on security paper using the VS Service Center letterhead stationery.